



**Ref.: OP1196D (right) / OP1196I (left)**



**SUPRAMALLEOLAR DYNAMIC ANKLE-FOOT ORTHOSES**

**Characteristics:**

- Dynamic orthoses with a supramalleolar design, made of flexible thermoplastic, comprising:
  - An internal part that surrounds the tibial and fibular malleoli and which comes into contact with the entire plantar surface. It embraces the entire perimeter of the midfoot reinforcing a proper configuration of the foot to the orthoses. An external part which serves to reinforce the alignment applied by the orthoses.
  - A closure system with velcro at the midfoot ensuring an adequate fitting. Some adhesive velcro pieces are included. These can be placed in the heel area between the two parts to avoid movement between the foot and the orthoses, favouring a comfortable and firm fitting.
  - It is very light, allows for plantar and dorsal flexion, realigns the entire foot and has a support for the internal longitudinal arch that facilitates the development of fundamental movement patterns for a more natural gait. To ensure greater comfort during use, pads for the toes and the dorsal midfoot are included.

**Indications:**

- Moderate or severe pronosupination.
- Instability in the subtalar joint, midfoot and forefoot as a result of low or high muscular tone.
- Severe flat foot.
- Lack of coordination or balance.

**Colour:**

- White.

SIZES	FOOT MEASUREMENT CM
1	12,4
2	13,2
3	14
4	15
5	16
6	17,5
7	19
8	21

It is recommended to choose a size that leaves approximately 1cm of space between the end of the orthoses and the foot.