

It stimulates and helps the step pattern.
It trains the user to maintain the balance of the trunk and pelvis.
It organizes a wider and structured vision of space, enhances the perception of distances and environment that surrounds the user.
It develops physical self-sufficiency and psychological independence by experimenting the possibility of acting independently without the help of other people.

Grillo compensates the user lack of balance thanks to its ergonomic support and wide base, it allows the user to concentrate his/her residual abilities on walking.

The wide range of components and supports offers maximum versatility. All accessories are interchangeable quickly and without the use of tools to suit different clinical pictures.

design and quality
MADE IN ITALY since 1980



Product registered
FDA USA

Product marked
CE



- **RESISTANT FRAME**, foldable for an easy transportation with **PATENTED HEIGHT ADJUSTMENT SYSTEM**.
- Padded, width and circumference adjustable **PELVIC SUPPORT**. Padded,

tilt, height, width, circumference and forward-backward adjustable **TRUNK SUPPORT**. They allow an easy access. They hold the user just like **HANDS** giving him/her a sense of security and stability allowing the user to concentrate on walking without any fear.

- **ERGONOMIC**, comfortable, adjustable and easily removable for washing **HARNESS**. It supports the user in case of necessity without generating painful pressures.
- Removable and height adjustable **HANDLEBAR**. It is used as a reference

and as an armrest.

- Exclusive **ORMESA DESIGN**, light and puncture proof **WHEELS** : 2 front swivelling Ø 17.5 cm and 2 rear fixed Ø 25 cm.
- **DIRECTIONAL LOCK** on front wheels. **ADJUSTABLE DRAG, ANTIREVERSE DEVICE, WHEELS STOP** on rear wheels.
- **SEMIELASTIC BUMPERS**.



Attention!
Read carefully the user and maintenance manual

FRONT DRIVE	POSTERIOR	USER DIMENSIONS:
<p>Supports circumference: size 5: cm 84-140 size 3-4: cm 74-120</p> <p>Weight kg. 23 MAX. LOAD kg. 110</p>	<p>Support circumference: size 5: cm 84-140 size 3-4: cm 74-120</p> <p>Weight kg. 22 MAX. LOAD kg. 110</p>	<p>USER HEIGHT (C): cm 150-190</p> <p>IMPORTANT: Make sure that measures A and B indicated in the figure do not exceed the suggested range indicated for the supports</p> <p>SUPPORTS DIMENSIONS: Pelvic support A (height cm) 63-117 Trunk support B (height cm) 81-143</p>

ORDER FORM

GRILLO Large SIZE L (size Five)

FRONT DRIVE gait trainer POSTERIOR gait trainer



ORMESA S.r.l. - Via A. Da Sangallo, 1
PACIANA - 06034 FOLIGNO (PG) ITALY
TEL. +39 0742 22927 FAX +39 0742 22637
info@ormesa.com - www.ormesa.com

DISTRIBUTOR/AGENT:

DATE OF ORDER:

DATE OF DELIVERY:

NAME OF THE ADDRESSEE/USER:

VERSIONS ▶

COMPONENTS ▼

925 and 924: YOU CAN CHOOSE
the SMALLER SIZE

865 Adjustable headrest		<input type="checkbox"/>
925 Trunk support		<input type="checkbox"/> M
924 Pelvic support		<input type="checkbox"/> M
930 Ergonomic harness		<input type="checkbox"/>
926 Handlebar		<input type="checkbox"/>
809 Arm supports		<input type="checkbox"/>
839 Ergonomic handles		<input type="checkbox"/>
815 Saddle		<input type="checkbox"/>
927 Removable grip handles		<input type="checkbox"/>
890 Leg divider with thigh straps		<input type="checkbox"/>
810 Weight bars kg. 3,4 x 2		<input type="checkbox"/>
928 Rear wheels for narrow base with pedals for parking brake (the frame width is reduced by 7.6 cm)		<input type="checkbox"/>
929 Pedals for parking brakes		<input type="checkbox"/>

<input type="checkbox"/> Front drive PT PELVIC-TRUNK-HANDLEBAR	<input type="checkbox"/> Front drive PA PELVIC-ARM SUPPORTS	<input type="checkbox"/> Posterior PT PELVIC-TRUNK SUPPORTS	<input type="checkbox"/> Posterior PA PELVIC-ARM SUPPORTS	<input type="checkbox"/> Posterior P PELVIC SUPPORT
INCLUDING: TRUNK SUPPORT, PELVIC SUPPORT, ERGONOMIC HARNESS, HANDLEBAR	INCLUDING: PELVIC SUPPORT, ERGONOMIC HARNESS, ARM SUPPORTS	INCLUDING: TRUNK SUPPORT, PELVIC SUPPORT, ERGONOMIC HARNESS, HANDLEBAR.	INCLUDING: PELVIC SUPPORT, ERGONOMIC HARNESS, ARM SUPPORTS.	INCLUDING: PELVIC SUPPORT, ERGONOMIC HARNESS, HANDLES.
NA	NA	<input type="checkbox"/>	<input type="checkbox"/> Available in combination with trunk support only	<input type="checkbox"/> Available in combination with trunk support only
•	<input type="checkbox"/>	•	<input type="checkbox"/>	<input type="checkbox"/>
•	•	•	•	•
•	•	•	•	•
•	<input type="checkbox"/> WITHOUT JOINT	NA	NA	NA
<input type="checkbox"/> WITHOUT JOINT	•	<input type="checkbox"/> WITH JOINT	•	<input type="checkbox"/> WITH JOINT
<input type="checkbox"/> WITHOUT JOINT	•	<input type="checkbox"/> WITH JOINT	•	<input type="checkbox"/> WITH JOINT
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
NA	NA	•	<input type="checkbox"/>	•
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

• = INCLUDED

☐ = ACCESSORY

NA = NOT AVAILABLE

Notes:



Congratulation for choosing **ORMESA** !