

Advance Bath Chair



<p>Base price: €</p> <p>All prices in EUR, valid from 01/01/2017</p>	<p>Order no.:</p> <p>-----</p>
--	--------------------------------

<input type="checkbox"/> Provide estimate only		<input type="checkbox"/> Order	
Billing address & Customer order number			
Delivery address			
Contact		Name	
		
		Date	Signature

Standard equipment

- Frame with angle-adjustable back and leg support made of stainless steel and aluminium and plastic parts.
- Removable, washable netting cover
- Hip belt of netting material
- Chest belt of netting material

Notice:

- Size 1: for children from 1 to 5 years /
body height from 75 to 105 cm
- Size 2: for children from 4 to 9 years /
body height from 95 to 125 cm
- Size 3: for children from 8 to 14 years /
body height from 115 to 145 cm
- Size 4: for children from 12 to 18 years /
body height from 135 to 165 cm

Leckey Advance Bath Chair standard model			
LAB/1/COMP-	<input type="checkbox"/>	Size 1, includes chest belt and hip belt	AC
LAB/2/COMP-	<input type="checkbox"/>	Size 2, includes chest belt and hip belt	AC
LAB/3/COMP-	<input type="checkbox"/>	Size 3, includes chest belt and hip belt	AC
LAB/4/COMP-	<input type="checkbox"/>	Size 4, includes chest belt and hip belt	AC

Colours				
AM	01	<input type="checkbox"/>	Aquamarine	-
RB	02	<input type="checkbox"/>	Royal blue	-
OR	04	<input type="checkbox"/>	Terracotta	-
YE	06	<input type="checkbox"/>	Sun yellow	-
PK	07	<input type="checkbox"/>	Pink	-

Accessories				
LAB/1/B-AM	415J80=SK010-01	<input type="checkbox"/>	Abduction belt for sizes 1 & 2, Aquamarine	AC
LAB/1/B-RB	415J80=SK010-02	<input type="checkbox"/>	Abduction belt for sizes 1 & 2, Royal blue	AC
LAB/1/B-OR	415J80=SK010-04	<input type="checkbox"/>	Abduction belt for sizes 1 & 2, Terracotta	AC
LAB/1/B-YE	415J80=SK010-06	<input type="checkbox"/>	Abduction belt for sizes 1 & 2, Sun yellow	AC
LAB/1/B-PK	415J80=SK010-07	<input type="checkbox"/>	Abduction belt for sizes 1 & 2, Pink	AC
LAB/2/B-AM	415J80=SK020-01	<input type="checkbox"/>	Abduction belt for sizes 3 & 4, Aquamarine	AC
LAB/2/B-RB	415J80=SK020-02	<input type="checkbox"/>	Abduction belt for sizes 3 & 4, Royal blue	AC
LAB/2/B-OR	415J80=SK020-04	<input type="checkbox"/>	Abduction belt for sizes 3 & 4, Terracotta	AC
LAB/2/B-YE	415J80=SK020-06	<input type="checkbox"/>	Abduction belt for sizes 3 & 4, Sun yellow	AC
LAB/2/B-PK	415J80=SK020-07	<input type="checkbox"/>	Abduction belt for sizes 3 & 4, Pink	AC
LAB/1/C-AM	415L80=SK010-01	<input type="checkbox"/>	Head support, Aquamarine	AC
LAB/1/C-RB	415L80=SK010-02	<input type="checkbox"/>	Head support, Royal blue	AC
LAB/1/C-OR	415L80=SK010-04	<input type="checkbox"/>	Head support, Terracotta	AC
LAB/1/C-YE	415L80=SK010-06	<input type="checkbox"/>	Head support, Sun yellow	AC
LAB/1/C-PK	415L80=SK010-07	<input type="checkbox"/>	Head support, Pink	AC
LSF/1	415A80=SK090	<input type="checkbox"/>	Shower frame for Leckey Advance Bath Chair	AC

Information



Shower frame

Technical Data

	Size 1	Size 2	Size 3	Size 4
Age of child	1 – 5 years	4 – 9 years	8 – 14 years	12 – 18 years
Body height	75 – 105 cm	95 – 125 cm	115 – 145 cm	135 – 165 cm
Seat depth	225 mm	275 mm	335 mm	400 mm
Headrest	150 mm	150 mm	215 mm	215 mm
Back support	350 mm	400 mm	460 mm	550 mm
Seat height	50 – 425 mm	50 – 425 mm	50 – 425 mm	50 – 425 mm
Overall width	430 mm	430 mm	430 mm	430 mm
Max. user weight	72 kg	72 kg	72 kg	72 kg

Shower trolley

	Size	Overall length (inches)
Width	600 mm	38
Length	900 – 1150 mm	43
Height	650 – 700 mm	52
	4	58