



# Leckey PAL

<b>Name</b>	
<b>Date:</b>	
<b>Email</b>	
<b>Contact Number</b>	
<b>Invoice Address</b>	
<b>Delivery Address</b>	

Quote       Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Price
<b>Seat</b>				
<input type="checkbox"/>	135-1600-03	PAL Size 1 Seat (including lap strap) - Blue		XXX
<input type="checkbox"/>	135-1600-01	PAL Size 1 Seat (including lap strap) - Green		XXX
<input type="checkbox"/>	135-2600-03	PAL Size 2 Seat (including lap strap) - Blue		XXX
<input type="checkbox"/>	135-2600-01	PAL Size 2 Seat (including lap strap) - Green		XXX
<input type="checkbox"/>	135-3600-03	PAL Size 3 Seat (including lap strap) - Blue		XXX
<input type="checkbox"/>	135-3600-01	PAL Size 3 Seat (including lap strap) - Green		XXX
<input type="checkbox"/>	135-4600-03	PAL Size 4 Seat (including lap strap) - Blue		XXX
<input type="checkbox"/>	135-4600-01	PAL Size 4 Seat (including lap strap) - Green		XXX
<b>Tray</b>				
<input type="checkbox"/>	135-1620	PAL Activity tray size 1		XXX
<input type="checkbox"/>	135-1656	PAL Tray Pad size 1		XXX
<input type="checkbox"/>	135-2620	PAL Activity tray size 2		XXX
<input type="checkbox"/>	135-2656	PAL Tray Pad size 2		XXX
<input type="checkbox"/>	135-3620	PAL Activity tray size 3		XXX
<input type="checkbox"/>	135-3656	PAL Tray Pad size 3		XXX
<input type="checkbox"/>	135-4620	PAL Activity tray size 4		XXX
<input type="checkbox"/>	135-4656	PAL Tray Pad size 4		XXX
<b>Side Pads</b>				
<input type="checkbox"/>	135-1635-03	PAL 25mm Side Pads size 1 (pair) - Blue		XXX
<input type="checkbox"/>	135-1635-01	PAL 25mm Side Pads size 1 (pair) - Green		XXX
<input type="checkbox"/>	135-1654-03	PAL 50mm Side Pads size 1 (pair) - Blue		XXX
<input type="checkbox"/>	135-1654-01	PAL 50mm Side Pads size 1 (pair) - Green		XXX
<input type="checkbox"/>	135-2635-03	PAL 25mm Side Pads size 2 (pair) - Blue		XXX
<input type="checkbox"/>	135-2635-01	PAL 25mm Side Pads size 2 (pair) - Green		XXX
<input type="checkbox"/>	135-2654-03	PAL 50mm Side Pads size 2 (pair) - Blue		XXX
<input type="checkbox"/>	135-2654-01	PAL 50mm Side Pads size 2 (pair) - Green		XXX
<input type="checkbox"/>	135-3635-03	PAL 25mm Side Pads size 3-4 (pair) - Blue		XXX
<input type="checkbox"/>	135-3635-01	PAL 25mm Side Pads size 3-4 (pair) - Green		XXX
<input type="checkbox"/>	135-3654-03	PAL 50mm Side Pads size 3-4 (pair) - Blue		XXX
<input type="checkbox"/>	135-3654-01	PAL 50mm Side Pads size 3-4 (pair) - Green		XXX
<b>Ramped Cushion</b>				
<input type="checkbox"/>	135-1651-03	PAL Ramped Seat Cushion size 1 - Blue		XXX
<input type="checkbox"/>	135-1651-01	PAL Ramped Seat Cushion size 1 - Green		XXX
<input type="checkbox"/>	135-2651-03	PAL Ramped Seat Cushion size 2 - Blue		XXX
<input type="checkbox"/>	135-2651-01	PAL Ramped Seat Cushion size 2 - Green		XXX
<input type="checkbox"/>	135-3651-03	PAL Ramped Seat Cushion size 3 - Blue		XXX
<input type="checkbox"/>	135-3651-01	PAL Ramped Seat Cushion size 3 - Green		XXX
<input type="checkbox"/>	135-4651-03	PAL Ramped Seat Cushion size 4 - Blue		XXX
<input type="checkbox"/>	135-4651-01	PAL Ramped Seat Cushion size 4 - Green		XXX
<b>Pommel</b>				
<input type="checkbox"/>	135-1602-03	PAL Pommel (size 1) - Blue		XXX
<input type="checkbox"/>	135-1602-01	PAL Pommel (size 1) - Green		XXX
<input type="checkbox"/>	135-3602-03	PAL Pommel (size 2 3 4) - Blue		XXX
<input type="checkbox"/>	135-3602-01	PAL Pommel (size 2 3 4) - Green		XXX
<b>Footrest</b>				
<input type="checkbox"/>	135-1603	PAL Footrest size 1		XXX
<input type="checkbox"/>	135-2603	PAL Footrest size 2		XXX
<input type="checkbox"/>	135-3603	PAL Footrest size 3		XXX
<input type="checkbox"/>	135-4603	PAL Footrest size 4		XXX
<b>Base</b>				
<input type="checkbox"/>	135-1610	PAL Ski with castors size 1 (pair)		XXX

<input type="checkbox"/>	135-1604	PAL Rocker / stabilisers size 1 (pair)		XXX
<input type="checkbox"/>	135-2610	PAL Ski with castors size 2 (pair)		XXX
<input type="checkbox"/>	135-2604	PAL Rocker / stabilisers size 2 (pair)		XXX
<input type="checkbox"/>	135-3610	PAL Ski with castors size 3 (pair)		XXX
<input type="checkbox"/>	135-3604	PAL Rocker / stabilisers size 3 (pair)		XXX
<input type="checkbox"/>	135-4610	PAL Ski with castors size 4 (pair)		XXX
<input type="checkbox"/>	135-4604	PAL Rocker / stabilisers size 4 (pair)		XXX

	Size 1	Size 2	Size 3	Size 4
Age	1-3	2-5	4-8	7-12
User Height	Min 700mm [27.6"] Max 950mm [37.4"]	Min 900mm [35.4"] Max 1050mm [41.3"]	Min 1000mm [39.4"] Max 1300mm [51.2"]	Min 1200mm [47"] Max 1600mm [63"]
Max. User Weight	30kg	35kg	40kg	50kg
Seat Depth	180mm-240mm [7"-9.4"]	220mm-290mm [8.6"-11.4"]	250mm-330mm [9.8"-12.9"]	300mm-400mm [11.8"-15.7"]
Seat Width	150mm-250mm [5.9"-9.8"]	170mm-270mm [6.6"-10.6"]	230mm-330mm [9"-12.9"]	280mm-380mm [11"-14.9"]
Seat Height	160mm-240mm [6.2"-9.4"]	185mm-285mm [7.2"-11.2"]	250mm-350mm [9.8"-13.7"]	300mm-400mm [11.8"-15.7"]