

tray and wheels included


DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FUNCTIONS

- **No edge tilting frame** in order to enable trunk rotation and side-arm motion. Smooth continuous adjustment, without snaps (through spring stopper adjuster). It's adjustable by knobs, without using wrenches, "tools free".
- **Birch wood table.** It is removable and adjustable in height, tilt and depth.
- **Trunk Support** adjustable in its circumference thanks to Velcro straps, in height, back and forward.
- **Pelvic Support** adjustable in its circumference thanks to Velcro straps, in height, back and forward.
- **Numbered bars** to find the right setting, determined by the therapist.
- **Knee-pads** adjustable in height, width, depth and in a rotatory movement.
- **Heelrests** adjustable in depth and width.
- **Birch wood footboard.** It gives a safe guide to the child. It is adjustable in height.
- **Twinned solid Wheels,** 75 mm. diam., all wheels with brakes. Thanks to the wheels the aid can be moved from one room to another having the child on it.



MARKED 
 FDA registered - USA


 **Attention!**
 Read carefully the user and maintenance manual

Please check off the interested products:



MEASURES

- DONDOLINO 1 DONDOLINO 2 DONDOLINO 3



User dimensions

USER HEIGHT (C)
 size 1: min cm 75 - max cm 100
 size 2: min cm 90 - max cm 125
 size 3: min cm 120 - max cm 150


IMPORTANT:
 make sure that measures F and G indicated in the figure do not exceed the suggested range indicated for the supports.

size	Max. Load Cap.	Weight	A	B	C	D	E	F	G	L	P
1	Kg. 35	Kg. 21	cm 63	cm 75-85	cm 45-79	cm 63	cm 12-26	cm 39-59	cm 50-76	cm 60	0°-19°
2	Kg. 45	Kg. 24	cm 77	cm 81-91	cm 52-90	cm 67	cm 12-36	cm 47-74	cm 56-90	cm 65	0°-14°
3	Kg. 55	Kg. 26	cm 83	cm 85-95	cm 62-107	cm 71	cm 14-48	cm 58-93	cm 70-110	cm 70	0°-11°

ACCESSORIES

944 Circumference reduction parts
 Thickness: 3 cm. each side

865 Multiadjustable headrest



Notes: _____

