

SCOOOT



3 in 1

The 3-in-1 mobility rider
for little explorers

Quick and easy ordering
www.fireflyfriends.com
Supporting Special Needs
Family Participation





① Crawl

Crawl lets kids move and explore on their tummy, helping to strengthen upper body and tummy muscles as they lift their heads and shoulders.



② Scoot

In Scoot function, kids can sit upright and move around using their feet – this encourages their trunk and pelvis to stabilise, freeing up their hands to handle objects, play and explore.

2-in-1 Scoot
Crawl and Scoot

3-in-1 Scoot
Crawl, Scoot and Ride



③ Ride

Ride lets kids self-propel and manoeuvre around the floor, building upper body strength and experiencing independent movement as they go.



NEW!
Advanced Backrest Cushion

This optional accessory offers kids additional support around the trunk, or if the laterals are removed, acts as a comfy cushioned backrest.

Sizing Guide	2-in-1 Scoot	3-in-1 Scoot
Approx. age	2-6 years	2-6 years
Maximum user weight	22kg/48.5lbs	22kg/48.5lbs
Seat base width (room for hips)	260mm/10"	260mm/10"
Seat base length (from back to front of seat)	400mm/15.5"	400mm/15.5"
Length from hip to end of footplate when seated	N/A When footplate is not extended: 520mm/20.5"	When footplate is fully extended: 630mm/24.5"
Backrest width	230mm/9"	230mm/9"
Backrest height	160mm/6"	160mm/6"
Wheel arch width (for fitting through doorways!)	N/A	400mm/15.5"
Height of seat base off floor	160mm/6"	160mm/6"

Advanced Backrest Cushion Accessory

Backrest width	Backrest height	Laterals minimum setting	Laterals maximum setting	Height of laterals
230mm/9"	160mm/6"	440mm/17"	680mm/26.5"	140mm/5.5"

Lap Belt
(not intended to provide support)

Minimum hip setting	Maximum hip setting
380mm/15"	560mm/22"