

# MINI STANDY

ENGLISH  
COD. 22081-03/07



ORMESA®

ORMESA®



ORMESA

Via A. Da Sangallo 1 Z.A.I. PACIANA  
06034 FOLIGNO - PG - ITALY  
TEL. +39 0742 22927 (r.a.) FAX +39 0742 22637  
E-mail: [info@ormesa.com](mailto:info@ormesa.com)  
<http://www.ormesa.com>

COMPANY  
WITH QUALITY MANAGEMENT  
SYSTEM CERTIFIED BY DNV  
=ISO 9001:2000=

A.D.M. - A.REHA  
ASSOCIATION OF MANUFACTURERS OF  
AIDS FOR MOTORIAL DISABILITIES



PRODUCT TESTED BY  
TÜV BERLIN



ADJUSTABLE SUPPORT with velcro strip.



Adjustable **KNEE PADS** for **FOOT FLEXION-EXTENSION**.



Adjustment of **ADDUCTION-ABDUCTION OF LEGS**.



Removable **WORK TABLE**. IT IS **ADJUSTABLE** in height and depth.



Adjustable **FRONT** and **BACK CHEST SUPPORT**.



Adjustable **FRONT** and **BACK PELVIC SUPPORT**.

**FRAME** with a **SMOOTH CONTOUR** to **ALLOW TRUNK ROTATION** and **SIDWAYS ARM MOVEMENT**.



If the child has **SUFFICIENT TRUNK CONTROL**, only the **PELVIC SUPPORT BAND** is necessary.

