


Rotation of the shoulder



Prono-supination of the elbow



Flexo-extension of the elbow



Deviation of the wrist



Technical features

Movements carried out on the shoulder:

- Elevation $5^{\circ} \div 180^{\circ}$
- Adduction-abduction $35^{\circ} \div 150^{\circ}$
- Internal-external rotation $90^{\circ} \div 0^{\circ} \div 90^{\circ}$

Movements carried out on the elbow:

- Flexo-extension $0^{\circ} \div 150^{\circ}$
- Prono-supination $90^{\circ} \div 0^{\circ} \div 90^{\circ}$

Movements carried out on the wrist:

- Flexo-extension $80^{\circ} \div 0^{\circ} \div 80^{\circ}$
- Ulnar and radial deviation $20^{\circ} \div 0^{\circ} \div 30^{\circ}$

Programmable functions:

- Speed in both directions: min. $1.5^{\circ}/\text{sec.}$ max. $5.5^{\circ}/\text{sec.}$
- Work time
- Automatic increase of the breadth of movement on both limits
- Pause in both limits of movement
- Warm-up cycles

Technical specifications:

- Power supply: $85 \div 260 \text{ V} \sim 50 \div 60 \text{ Hz.}$
- Electric safety: Class 1 B Standard EN 60601-1
- Electromagnetic compatibility: Group 1 Class B Standard EN 60601-1-2
- Classification according to EEC Directive 93/42: Class IIa
- Weight: 77 Kg.

The instrument is compliant with the applicable Basic Requirements for safety, Health and protection established by EEC Directive 93/42.



Flexo-extension of the wrist



THE ACCESSORIES


Fisiotek HP2 in its base model allows to rehabilitate the shoulder.

The complete set of accessories (FHP229) includes the accessories to rehabilitate the elbow and the wrist.

The set of accessories Elbow (FHP230) includes the accessories to rehabilitate the elbow.

The set of accessories Wrist (FHP233) includes the accessories to rehabilitate the wrist.


Adjustable rod for arm support
Cod. FHP221




Fixed elbow support
Cod. FHP227



Handgrip for prono-supination of the elbow
Cod. FHP225




Trolley for accessories
Cod. FHP228



Support for ulnar and radial deviation of the wrist
Cod. FHP226




Handgrip for flexion-extension of the elbow
Cod. FHP224



Fixt wrist support
Cod. FHP223



Handgrip for flexion-extension of the wrist
Cod. FHP222



The seat accessory is also available
Cod. FHP236

MEDIOTRONICS (Pty) Ltd

545 Nupen Crescent,

Halfway House, Midrand

P.O. Box 7709, Halfway House, 1685

Tel: 011-315-2619

Fax: 011-315-2242

Web Address: www.mediotronics.co.za

Rimec S.r.l.

Loc. Braine, 57/a

I-40040 Rioveggio (Bologna) Italy

Tel. ++39-051-67.77.798

Fax ++ 39-051-67.77.357

E-mail: rimec@rimec.it

www.rimec.it