

Smart TENS

With Han Stimulation*

This powerful state-of-the art digital TENS machine is smooth, sophisticated and simple to use. It has 4 pre-set modes for ease of use. Ideal for treating numerous pain conditions.

"I bought the Smart TENS as I wanted something stylish and simple. It is absolutely brilliant and is so smooth and I have no problems using it at all."

Sian Thomas



CE
0120


Key Benefits

- Dual channel – treats 2 pain areas
- Large digital display
- 4 preset modes; continuous, burst, modulation and mixed frequency (Han Stimulation)
- Pulse width adjustment
- Timer function
- Ergonomically designed
- Neck cord
- Rotating belt clip
- Size: 11.5 x 5.5 x 3cm
- Weight: 100g



What's included

- Pack of 4 electrodes
- Instruction manual
- 2 leadwires
- Carrying pouch
- Neck cord for mobility
- 2 x AA batteries



**Professor Han suggests that switching between low and high-frequency stimulation for three seconds each i.e. 2/100 stimulation produces a simultaneous activation of the of the enkephalin and dynorphin systems which creates a more potent analgesic effect to combat pain.*

Published in the January issue of 2003 of "Trends in Neuroscience." (26:17-22)



Order on line: www.bodyclock.co.uk or telephone 0500 22 00 61