

# JELOSIL SA

## INFRA 750



THERAPEUTIC LAMP WITH INFRARED RAYS FOR  
MEDICAL USE ACCORDING TO EUROPEAN STANDARD  
EN 60601 - 1 CLASS I TYPE B, AND EEC DIRECTIVE  
93/42 "MEDICAL DEVICE" CLASS IIa.

Cod.OA1B100130

The device is equipped with 3 IR lamps of 250W 230V each one, placed on a base with wheels. It's possible to regulate the inclination of the lamps .

This device can be programmed for a functioning until 15 minutes. When the programmed time is over, there is an automatic switching off.

Each IR lamp has a switch for independent switching on/off

#### APPLICATIONS:

Use in orthopedic field: treatment of arthritis, lumbago, rheumatism, sprains, and generally where a penetrating heat source is required.

Use in physiotherapeutic field: massage preparation and body heating for a better absorption of creams by the skin.

Heating: in cold rooms it's possible to obtain an immediate heating for the irradiated person.

FOR THE USE OF THIS DEVICE IT'S NECESSARY TO FOLLOW CAREFULLY

THE DOCTOR'S INSTRUCTIONS ABOUT EXPOSITION TIMES,  
APPLICATION DISTANCE AND FREQUENCY.

#### TECHNICAL INFORMATION:

Voltage:	230 V; 50 Hz
Power:	750W
Fuses:	6,3A
Lamp dimensions:	cm. 55x15x21
Regulating height:	from cm. 82 to cm.152
Weight:	Kg. 14
<u>ENCLOSED ACCESSORIES</u>	
- N. 2 pair of protective goggles	cod 0R0A100500
<u>SPARE PARTS</u>	
- 250W IR bulb - bulb - socket type E/27	cod 0R1A101640
- Fuses type 20x5 F 6,3A 250V	

**CE** 0051

Rev. E05 of 13/04/2010