

corner sitter



LECKEY®



This portable **Corner Sitter** unit with its optional table is an invaluable aid for both carers and users.


Corner Sitter Features and Accessories

Standard Features

- a Adjustable sides for increased lateral support
- b Angle adjustable backrest for both prone and recline
- c Adjustable depthommel
- d 45 degree hip belt
- e Integral handle for ease of carrying

Optional Accessories

- f Multi adjustable Contoured High back
- g Height and angle adjustable activity tray



Fabrics

The Corner Sitter comes in 2 fabrics.
Size 1 comes as standard in Farmyard fabric to suit the nursery environment. Size 2 and 3 in Jazz fabric.



Prone and recline angled backrest



Product Sizing

		Size 1		Size 2		Size 3	
Age (approx)							
min		1.5		5		8	
max		6		9		14	
		Size 1		Size 2		Size 3	
		mm	inches	mm	inches	mm	inches
Overall User Height		1165	45.9	1330	52.4	1630	64.2
Seat Base							
length		335	14	460	18.1	550	21.6
width		380	15	425	16.7	475	18.7
Height in Folded Position		180	7	180	7	180	7
Weight of Product		9kg		11kg		13kg	
Backrest Height		300	11.8	400	15.7	500	19.7
High Back Height		120	4.7	140	5.5	160	6.3
High Back Adjustment		170	4	170	5.3	170	6.7
Pommel Adjustment		50	2.7	75	2.7	75	2.7
Max User Weight		30kg		40kg		50kg	
Degree of Angular Adjustment on Backrest							
prone		15		15		15	
recline		20		20		20	
Tray Sizes		640	25	640	25	640	25
		x450	x17.7	x450	x17.7	x450	x17.7
Tray Height (from ground)							
min		220	8.8	220	8.8	220	8.8
max		320	12.8	320	12.8	320	12.8



LS35-02

Distributed by



A member of the **LECKEY**® Group

_James Leckey Design Ltd
 _Kilwee Business Park
 _Dunmurry
 _BT17 OHD
 _Northern Ireland
 _United Kingdom

_T: 028 9060 2277

_freephone: _F: 028 9060 0795

_United Kingdom 0800 318265 _E: info@leckey.com

_Ireland 1800 626020 _W: www.leckey.com