
LECKEY°



LECKEY°

Stronger together.

We work with individuals, therapists and carers to design products with both a clinical and an emotional focus. Using the latest research and clinical understanding, we create practical solutions which are easily integrated into family life, because for us, life is about going, enjoying, participating and doing.

24 hour
postural care for
babies, kids & adults.
Sleeping, Sitting,
Standing, Walking,
Moving, Bathing,
Toileting.



Advance Bath Chair

The Advance Bath Chair is available in 4 sizes for infants through to young adults. The product provides the user with appropriate levels of postural support ensuring they are safe and secure at all times during bathing.

The simple and easy to use adjustments mean that the product can be set up in a variety of positions in seconds depending on the user's postural requirements.

The height adjustable legs mean that the product can be set up at the appropriate height for the parent or carer to reduce the risk of back strain or injury.

Bath Support

The Bath Support, available in 2 sizes for children aged 1-9 years of age. complements the Advance Bath Chair, providing support for the less involved child.

The bath support easily folds flat for storage & transportation.

Shower Trolley

With the addition of the adjustable Shower Trolley, the Advance Bath Chair can be converted in seconds to a versatile shower chair.

Manufactured from medical grade stainless steel, the robust frame and high quality lockable castors ensure the trolley can be easily and safely manoeuvred over wet and slippery floors.











inches

12 - 18

mm

1350 1650

LAB/4



inches

Product Sizing

Size	1		2		3
Age	1 - 5		4 - 9		8 - 14
	mm	inches	mm	inches	mm
User Height					
Min Max	750 1050	30 41	950 1250	37 49	1150 1450
Advance Bath Chair	LAB/1		LAB/2		LAB/3
Head support Back support Seat depth Calf support Seat Height 50-425mm / 2-17" Max Weight 72kg / 159lbs	150 350 225 260	6 14 9 10	150 400 275 260	6 16 11 10	215 460 335 340
Total Width 430mm / 17"					
Shower Trolley	LSF/1				
Width Length:	600	24			
Min Max	900 1150	36 45			
Seat Height from floor 650-700mm / 26"- 28"					
Bath Support	LBS/1		LBS/2		
Head support Back support Seat depth	150 600 300	6 23.6 12	150 700 400	6 27.5 15.7	
Seat Height 0-150mm / 0-6" Max Weight 62kg / 137lbs Total Width 400mm / 16"					

Colour Options

The fabric, available in seven bright colours, is rinse clean, breathable, soft to the touch and can be machine washed at 60°C.

















Accessories



LTS/#/C Depth adjustable tray



LAB/1/A Chest belt with laterals OTS/#/A Mounting bracket



Toilet Seat

The Toilet Seat comes in three sizes to suit the needs of a wide age range of children and adults.

It can be used as a: Freestanding Commode; Shower Chair; Over Toilet Seat.

The Toilet Seat is extremely versatile offering both toileting and shower options. As with our Bathing Range, we set out to develop products using the highest quality materials, delivering designs which are hygienic, safe and durable.

The padded backrest can be set to three different positions to reduce or increase the seat depth to ensure the user has a stable base of support.

The Toilet Seat is made from high performance wipe clean polymers which are pre-treated with a non-toxic, anti-fungal and anti-bacterial agent.

Over Toilet Seat

The Toilet Seat can be easily mounted onto any standard toilet when fitted with the over-toilet mounting bracket.

The seat can be simply removed from the standard toilet, as and when required.





Over toilet seat



Mounting bracket

Product Sizing

Size	1		2		3				
Age	1 - 5		4 - 9		8 - 14				
Codes	mm	inches	mm	inches	mm	inches			
Toilet Seat Over Toilet Seat	LTS/1 OTS/1		LTS/2 OTS/2		LTS/3 OTS/3				
Backrest	400	16	500	20	600	24			
Seat Depth Min Max Seat width Seat to tray	180 240 220 200	7 9.5 8.5 8	220 330 290 250	8.5 13 11.5	240 400 370 300	9.5 16 14.5 12			
Seat to Footplate Min Max	160 400	6 16	160 400	6 16	160 400	6 16			
Backrest Angle	Prone / Re	Prone / Recline 10°							
Seat to footplate	160-400m	160-400mm / 6-16"							
Tray	W 500mr	W 500mm / 20" D 330mm / 13"							
Castor Frame	W 520mr	W 520mm / 20.5" H 450mm / 17.5" L 625mm / 25"							

Easy Seat Potty Trainer



Accessories



ht adjustable tray Seat pad for Easy Seat



Potty insert for Fasy S



asy Seat Footpl



With sandal straps

Product Sizing

Size	1		2		3		4		
Age	1-3		2 - 5	2 - 5		3 - 8		6 - 11	
Codes Fixed Legs Castor Legs	LPT/1/FL LPT/1/CL		LPT/2/FL LPT/2/CL			LPT/3/FL LPT/3/CL		LPT/4/FL LPT/4/CL	
	mm	inches	mm	inches	mm	inches	mm	inches	
Easy Seat Potty Trainer Length Width Weight of Product	400 350 12kg	16 14	480 370 15kg	19 14	570 450 20kg	22 18	620 490 22kg	24 19	
User Standing Height min max	760 1000	30 40	930 1150	37 46	1000 1300	40 52	1150 1450	46 58	
Max User Weight	30kg	66lbs	30kg	66lbs	40kg	88lbs	40kg	88bs	
Backrest Height (sitting shoulder height) min max	225 285	9 11	265 340	10 14	325 425	13 17	400 500	16 20	
Seat Base Depth (backrest to back of knee) min max	175 250	7 10	225 300	9 12	250 350	10 14	300 400	12 16	
Seat Base Width (hip width) min max	125 225	5 9	150 250	6 10	225 325	9 13	275 375	11 15	
Calf Length (sole of foot to back of knee) min max	150 230	6 9	220 300	9 12	250 350	10 14	300 420	12 17	
Tray Height (seat to top of tray) min max	135 225	5 9	135 225	5 9	180 325	7 13	180 325	7 13	
Armrest Height (seat base to elbow) max min	110 150	4 6	120 160	5 7	160 210	6 8	190 240	7 10	
Degree of angular adjustment on backrest Prone Recline	5° 10°		5° 10°		5° 10°		5° 10°		
Castor size	75	3	75	3	75	3	75	3	
Tray sizes	450 x420	18 x16	450 x420	18 x16	590 x520	23 x20	590 x520	23 x20	

Please refer to **www.leckey.com** for our most up to date product codes

Leckey Servicing

Our FREE service and support includes: Free product training Free product set up Free product assessment Free product re-assessment Free repair within 2 year warranty For service policies on all products outside warranty, please contact Leckey's Customer Service department. James Leckey Design Ltd as manufacturer with sole responsibility declares that all products conform to 93/42/EEC guidelines and EN12182 technical aids for disabled persons general requirements and test methods. Order forms and spare parts lists to extend the service life of the product and allow reissue are available on request or online at www.leckey.com.

LECKEY°

Established in 1983, Leckey is a globally recognised pioneer in the research and development of products that help adults and children with disabilities to go, do, enjoy and participate in everyday activities throughout the day and night.

We take a highly clinical approach to product design and development. Through in-depth clinical research studies with leading universities, and extensive trials with occupational therapists, physiotherapists, users and their families, we continue to develop posturally supportive, family friendly products for all day care, at every stage of life.

Through early intervention, childhood and adulthood Leckey's experienced team of designers, therapists and bioengineers work together to develop products that meet the clinical needs of the healthcare professionals and the social needs of the user.

To achieve this, we work with the healthcare professionals, the individuals and carers who use our products everyday. With their help, we create the dependable, durable, proven and high performance products that we are known for worldwide.









































Leckey 19 Ballinderry Road Lisburn **BT28 2SA** Northern Ireland United Kingdom

T: 028 9260 0750 F: 028 9260 0799 E: info@leckey.com W: www.leckey.com





24 hour postural care for babies, kids & adults. Sleeping, Sitting, Standing, Walking, Moving, Bathing, Toileting.