

# madita-fun.

# Always ready for some fun.

**GMFCS**  
Levels III-V

kg  
max. 50



(-)5°–25°

(-)5°–37°

15–42 cm

25–62 cm

16–38 cm



## I am your integrative therapy chair.

Whether at the breakfast table or in the play corner; I will put you on a level with your friends. I will move you up and down. It is easy due to the support from my gas pressure springs. Use my lever to tip me forwards and backwards.

I can be adapted to your body and your tastes: I am available in lots of great colours. Did you know that my backrest can be removed? This trains your back musculature and helps you to sit up straight.



**Push handle**

for easy transfer.

**Back unit**

optionally removable, can be adjusted in angle and height.

**Foot support**

for good positioning of the feet.

**Castors**

available in different sizes, facilitate transfer.

**Headrest**

can be adjusted in height, depth and angle.

**Lateral pads**

allow a stable and upright sitting position.


**Therapy table**

offers a lot of space for activities.

**Seat surface**

is individually adaptable.

**Foot pedal**

for height adjustment.

**Base frame**

compact and sturdily built.



**In just one single step, my back unit can simply be removed. This gives you more freedom to move.**


**Frame colours**

● Anthracite metallic

● Telegrey 4





Data	Size 0	Size 1	Size 1b	Size 2
Seat depth	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
Back height with short back	35 - 43 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
Back height with long back	42 - 52 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
Seat height	25 - 62 cm	25 - 62 cm	25 - 62 cm	25 - 62 cm
Seat height with electric height adjustment	32 - 62 cm	32 - 62 cm	32 - 62 cm	32 - 62 cm
Seat height with hydraulic height adjustment	27 - 63 cm	27 - 63 cm	27 - 63 cm	27 - 63 cm
Lower leg length	17 - 24 cm	22 - 32 cm	27 - 38 cm	30 - 48 cm
Total width	54 cm	54 cm	64 cm	64 cm
Length of base frame	61 cm	61 cm	75 cm	75 cm
Total length with combination sub-frame	76 cm	76 cm	86 cm	86 cm
Seat angle	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°
Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
max. load	30 kg	30 kg	45 kg	50 kg
Weight	19 kg	20 kg	22 kg	24 kg
HMV No.	26.11.05.0022			



You can find more information on fabrics & embroidery on pages 188 - 191.

You can find more information on accessories & harnesses in the price list.

