

GMFCS  
Levels I-III

kg  
max. 60



(-)3° -17.5°

(-)5° -15°

16-40 cm

17-40 cm

### I am your seating aid for everyday use.

I promote active sitting and thus support your development. The great thing is whatever size you are, I will grow with you! You can simply adjust my seat and footplate up or down or to the rear so that you can sit comfortably at the table with your big brother. I have been manufactured with loving attention to detail using wood. I suit the other furniture in your home perfectly! With my many accessory parts, I can easily be adapted even more individually to your requirements.





**You decide how I should look.  
Combine my frame in natural  
beech or white stain  
with a  
selection of great coloured  
covers!**



**Headrest**

supports the head control.

**Backrest**

can be adjusted in height and angle.

**Armrests**

can be adjusted in height.

**Seat plate**

can be adjusted in height, depth and angle.

**Base foot**

with plastic sliders.

**Removable pad cover**

made from breathable material.

**Lateral pads**

allow a stable and upright sitting position.

**Therapy table**

offers a lot of space for activities.

**Thigh guide**

offers necessary support.

**Wooden frame**

with integrated aluminium profiles.

**Footplate**

can be adjusted in height, depth and angle.

**Fast height adjustment**

via the aluminium profiles integrated in the frame.

**Fast angular adjustment**

via pull strap





With a great  
range of  
accessories!



Data	Size 0 Nursery	Size 0	Size 1 Nursery	Size 1	Size 2
Seat width	17 – 28 cm	17 – 28 cm	20 – 35 cm	20 – 35 cm	30 – 40 cm
Seat depth	16 – 26 cm	16 – 26 cm	25 – 35 cm	25 – 35 cm	30 – 40 cm
Max. seat height	45 cm	65 cm	45 cm	65 cm	65 cm
Back height	25 cm	25 cm	29 cm	29 cm	34 cm
Back angle	(-)3° – 17.5°				
Seat angle	(-)5°, 0°, 7.5°, 15°				
Foot angle	(-)5°, 0°, 7.5°, 15°				
Foot support height min.	13 cm* / 18 cm**/ 21 cm***				
Foot support height min. (with fast angle adjustment)	8 cm* / 13 cm**/ 16 cm***				
max. load	30 kg	30 kg	45 kg	45 kg	60 kg
Weight	9.2 kg	9.7 kg	9.6 kg	10.2 kg	12 kg
HMV No.	26.11.05.0034				

Adjustable lower leg length max. = selected/required seat height – foot support height min. \* Standard, \*\* With 50 mm double castor, \*\*\* With 75 mm double castor



You can find more information on fabrics & embroidery  
on pages 188 – 191.

You can find more information on accessories & harnesses  
in the price list.

