Classroom Seat



Posture + Attention = Learning

Designed by the Leckey team of experienced therapists and engineers, Pal is an attractive, robust and easy to use classroom and nursery seating system for kids of 1-12 years of age, with mild to moderate postural needs. It is designed to give improved stability, which helps reduce fatigue and allows the child longer periods of concentration and fine motor activity. Available in 4 sizes and 4 colours it will help kids up to 12 years of age to participate, learn and give their full attention in the classroom.







Created for self assessment

The simplicity of Pal enables the therapist to assess the child's needs, instantly order their Pal Classroom Seat and quickly and easily set up the product when it arrives. Pal really is the simple, fast way to help kids join in with learning.



Available in 4 sizes and 4 colours 12344




Order using 4 easy steps...

1. Measure

Pal is a self assessment product where you assess the size required by

- 3 simple measurements:
- 1. Seat depth
- 2. Seat width
- 3. Seat height



Pal is available in 4 sizes and comes in 4 fab colours: green, orange, blue and pink.

2. Order

Email info@leckey.com or just call our Customer Service Department **UK** 0800 318 265 ROI 1800 626 020 for an instant quote indicating:

1. Chair size and colour.

2. Any leg accessory: Castor, Rocker, Stabiliser or High Seat.

3. Any other accessory items: Footrest, Tray, Table, Pommel and the Side Pad option you require.

3. Delivery

Our fast track service will deliver your chair within 14 days. The chair will be delivered assembled and ready for use.

4. Set-up

Set the leg height, armrest height and back angle to suit the child. Easy!!



Pal Seat with rocker leg accessory

Pal Seat with high leg accessory (Size 1 only)

Pal Sizing Chart

	Age kange	Stature Range	User Weight	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Seat Weight	Table Height
Size 1	1 - 3 years	700 - 950mm	Max - 30kg	180 - 240mm	150 - 250mm	160 - 240mm	120-150mm	265mm	4.3kg	400mm
Size 2	2 - 5 years	900 - 1050mm	Max - 35kg	220 - 290mm	170 - 270mm	185 - 285mm	140-170mm	315mm	6kg	460mm
Size 3	4 - 8 years	1000 - 1300mm	Max - 40kg	250 - 330mm	230 - 330mm	250 - 350mm	180-230mm	380mm	8.2kg	530mm
Size 4	7-12 years	1200 - 1600mm	Max - 50kg	300 - 400mm	280 - 380mm	300 - 400mm	190-240mm	450mm	10kg	n/a

The Pal Seat comes as standard with depth adjustable seat base, angle adjustable back rest, 2 point pelvic strap, height adjustable legs, height adjustable armrests and wipe clean back rest and seat base cushions.

Pal Product Codes

X = Size 1, 2, 3 or 4 AA = Colour (01 Green, 02 Orange, 03 Blue, 04 Pink)



135-X600-AA



Ski with Castors (Pair) 135-X610



Rocker (Pair) 135-X604



Stabilisers (Pair) 135-X604



High Seat Frame (Size 1 only) 135-1630



Footrest 135-X603




25mm Side Pad (Pair) 135-X635-AA **50mm Side Pad (Pair)** 135-X654-AA



Depth Adjustable Activity Tray 135-X620 Tray Pad (Not shown) 135-X656



Pommel Size 1 - 135-1602-AA Size 2 & 3 - 135-3602-AA



Ramped Seat Cushion 135-X651-AA



4 Person Table 1300x500mm 135-X701-AA







Leckey 19 Ballinderry Road Lisburn BT25 2SA Northern Ireland United Kingdom T: 02892 600750 F: 02892 600799 E: info@leckey.com W: www.leckey.com

LS210-03