



# GoTo Seat

One Little Seat, Countless Possibilities

A lightweight, portable and stylish postural support seat that enhances children's lives through participation.



**Product Features**

- Suitable for attachment to most standard chairs in the home, restaurants and classroom by using the adjustable attachment straps. Can also be attached to shopping trolleys.
- Adjustable headrest and laterals that adjust for growth.
- Standard or Advanced Headrest available.
- Pressure relief cushion available.
- Optional floorsitter accessory has 26 recline positions, from upright (106°) to right back (140°).



Available in 3 fabric colours and a vinyl version

Technical Details		
	Size 1	Size 2
Age	1-3	3-8
Seat Base to Top of Headrest	380mm-540mm (14"-21")	510mm-680mm (20"-27")
Lateral Support Width	190mm-260mm (7.5"-10")	210mm-300mm (8"-12")
Seat Depth	175mm (7")	244mm (9.5")
Max. User Weight	15kg (33lbs)*	30kg (66lbs)**

\*If a child is above 13kg/28.5lbs GoTo Seat cannot be used in a shopping trolley.  
 \*\*Size 2 GoTo Seats are too large for use in shopping trolleys.

