

firefly **Scooot**
Stylish Wheels for Little Explorers

A 4-in-1 sit and ride which gives children the freedom to explore and have fun in style.



Product Features

- Scoot has four configurations: Scoot, Ride, Pooosh and Crawl. Different configurations for different abilities and different ways to play.
- For added support in Scoot, the Advanced Backrest Cushion helps. It's an optional extra that gives support to the child's trunk. The laterals move up and down to give support where it's needed. Or remove them and just have a comfy cushioned backrest. The Advanced Backrest Cushion is easy to attach, wash, and can be bought separately.
- One size caters for children aged 2-6- the footrest adjusts to suit the child's leg-length as they grow.



Four Configurations

- **Crawl** allows kids to move and explore on their tummy, helping to strengthen their upper body and tummy muscles as they lift their heads and shoulders.
- **Scoot** lets kids sit upright to move around using their feet. Supported by a lap belt and backrest, the child's trunk and pelvis are encouraged to stabilise, freeing up their hands to lift and handle objects, play with family, and explore their surroundings.
- **Ride** lets kids propel themselves around the floor by moving the wheels with their hands. Having fun in this function can build upper body strength, whilst exploration and participation can improve spatial awareness as well as cognitive, sensory and communication abilities.
- **Pooosh** lets somebody else do the hard work. If the child isn't ready to propel themselves yet (or they're just tired), attach the Pooosh handle. It helps get little ones used to the Scoot. And it means parents and therapists can help them follow the fun at home, school or therapy sessions.

| Technical Details | |
|-------------------------|---------------------------|
| Size | 4-in-1 |
| Age | 2-6 |
| Max. User Weight | 22kg (48.5lbs) |
| Seat Depth | 400mm (15.5") |
| Seat Width | 230mm (9") |
| Hip to end of footplate | 520mm-630mm (20.5"-24.5") |
| Backrest Width | 260mm (10.5") |
| Backrest Height | 300mm (12") |
| Wheel Arch Width | 230mm (9") |
| Pooosh Pole | 800mm-890mm |

