

firefly **Splashy**

Water Adventures for Little Splashers

A lightweight, portable multipurpose bath seat that can be used for both bathing and other fun water play activities.



Product Features

- Splashy is designed to sit low in the bath to allow the child to benefit from being immersed in the water.
- Great aesthetics via bright colour options and a modern design.
- Splashy adapts to the child’s postural support needs, with the option of a 5-point harness or the 3-point pelvic support. There are four detachable bumpers to add extra support wherever it’s needed. The harness and foam supports attach easily to the seat base using toggles.
- Warm-touch foam feels soft whether it’s wet or dry.
- Splashy’s base unit lets you switch easily between 26 recline positions, from upright (106 degrees) to full recline (140 degrees).
- Lightweight, detachable parts make it easy to clean.



Available in 3 colours



Technical Details	
Age	1-8
Backrest Height	610mm (24")
Backrest Width	385mm (15")
Backrest Depth	300mm (12")
Seat Width	350mm (14")
Floorsitter Width	420mm (16.5")
Floorsitter Length	630mm (24")
Max Weight	30kg (66lbs)



firefly **SplashyBig**

Water Adventures for Big Splashers

With all of the features and benefits of Splashy, the SplashyBig opens up many more opportunities for 8-16 year olds who want to participate in fun activities with friends and family.

Product Features

- Suitable for 8-16 years with a max user weight of 55KG/122lbs.
- SplashyBig is lightweight and super easy to maintain with its warm touch and quick drying covers. Cleaning it down after use is a breeze.
- SplashyBig's floorsitter has 26 recline options plus a new stabiliser feature for use outside of the bath.



Available in 2 colour combinations

Technical Details

Seat Backrest Height	720mm [28.5"]
Seat Backrest Width	450mm [17.5"]
Seat Base Depth	400mm [15.5"]
Seat Base Width	450mm [17.5"]
Floorsitter Width	420mm [16.5"]
Floorsitter Width with Stabiliser	665mm [26"]
Floorsitter Length	795mm [31"]
Maximum Weight	55kg [121lbs]
Age Range	8-16 years old

