



Upsee

Little Steps to an Unexplored World

A mobility harness designed to allow children with special needs to experience upright movement like anyone else.



Product Features

- Facilitates gradual weight-bearing and gives children with motor impairment the closeness and confidence needed to stand and move with the help of an adult.
- Creates opportunity for social and emotional development as users can interact with peers at eye-level.
- Comprised of three parts - child harness, adult belt and double sandal.
- Durable, machine washable and breathable fabrics.
- Two colour options available for each size - please visit our website for further details.



Available in 4 fabric colours

Technical Details				
	Extra Small	Small	Medium	Large
Age	1-2	2-4	4-6	6-8
Max. User Weight	15kg (33lbs)	15kg (33lbs)	20kg (44lbs)	25kg (55lbs)
Shoulders to Hip	280mm-320mm (11"-13.5")	310mm-360mm (12"-14")	350mm-390mm (14"-15.5")	380mm-420mm (15"-16.5")
Chest Circumference	500mm-570mm (19"-22")	520mm-600mm (20"-24")	570mm-650mm (22"-26")	630mm-710mm (25"-28")
Lower Hip Circumference	500mm-570mm (19"-22")	520mm-600mm (20"-24")	580mm-670mm (23"-26")	640mm-730mm (25"-29")

