


# Travel SIT

ORTHOPAEDIC  
SEAT

NEW  
2019

LIGHT  
MOBILE  
COMPACT



WEIGHT : 3,7 KG

RANGE : WORLDWIDE

PREPARATION TIME : 30 SECONDS

HAPPINESS & INTEGRATION OF THE CHILD : 100%

RANGE OF COLOURS for Travel SIT





Travel SIT is light and mobile system for child stabilisation in sitting position, lying position and intermediate positions. The seat fits to the comfortable, dedicated for it backpack, which from now on will be your number one fellow traveller.

Backrest consists of four smoothly adjustable trunk supports, which together with the two safety belts and anatomical headrest has been made to support your child. Seat is equipped with six supports which are enabling stabilising child and adjusting the width of the seat. In the set you can also find adaptable abductor – which helps with the lower limbs abduction – and pelvic belts. Additional equipment is a lounge frame, which is enabling for the Travel SIT to be used on the ground.

This and more is available with the Travel SIT!


Start your journey, and be sure that Travel SIT will take care for proper sitting position of your child.



Integration with siblings and peers.



Share dinner at the restaurant.




Fun on the swing at the playground.



Fun next to the garden table.



Stay at the grandparents place.



Lack of the access to the special equipment during the holidays? Thanks to the TravelSIT – problem is solved!

## TECHNICAL DATA

	Size S	Size M
Patient's height	70÷105 cm	95÷125 cm
Backrest height	52 cm	62 cm
Seat's depth	26 cm	31 cm
Seat's width	17÷24 cm	22÷35 cm
Side supports width range	16÷22 cm	23÷30 cm
Side supports height range	8÷26 cm	13÷35 cm