



Technical information	Size 1	Size 2	Size 2+ *	Size 3	Size 4
User height	60 - 95 cm	90 - 130 cm	90 - 130 cm	120 - 155 cm	130 - 180 cm
User max weight	30 kg	40 kg	40 kg	55 kg	80 kg
Manual tilt	Yes	Yes	Yes	Yes	No
Power tilt	Yes	Yes	Yes	Yes	Yes
Footprint (length x width)	75 x 58 cm	75 x 58 cm	100 x 70 cm	100 x 70 cm	108 x 75 cm
* ELLA size 2+ , larger base allows maximal abduction					

MANUFACTURER

Timoteos Oy
Palorannantie 40
28660 Pori, FINLAND

VAT: FI10408720
email: timoteos@timoteos.fi

INTERNATIONAL SALES


Mikko Pirkola
tel. +358 50 385 6203
email: sales@timoteos.fi



Knee flexion



Tilt and rotation function available for pelvic and thoracic pads



Wide range of modules



Extra abduction (4-in-1)

ELLA STANDER

- Top class modularity; all modules fit every frame size
- Supports with stepless and independently bilateral adjustments
- Designed and manufactured in Finland
- Intuitive and simple adjustments make adaptation incredibly easy.
- 90° tilt with gas spring or electric actuator
- Upright, prone, supine and abducted standing, all with one product

MANUFACTURER

Timoteos Oy
Palorannantie 40
28660 Pori, FINLAND

VAT: FI10408720
email: timoteos@timoteos.fi

INTERNATIONAL SALES

Mikko Pirkola
tel. +358 50 385 6203
email: sales@timoteos.fi