


Hip flexion



Versatile adjustments for legs and feet



Anatomically correct hinge point positioning



Tilt and rotation adjustability available for lumbar and back support

## ODE SEATING SYSTEM

- Optimal, individual positioning
- All adjustments for size and functionality placed outside the seat. No need for excess client transfers, removing padding or replacing parts while making adjustments
- Hip flexion adjustment mechanism as standard
- Left and right side accurately and separately adjustable. Stepless, bilateral adjustments
- Optional dynamic back rest, with easy adjustment and locking
- Versatile adjustments for legs and feet as standard

Technical data	Mini	Size 1	Size 2
User max weight	25 kg	40 kg	55 kg
Seat depth	15 - 24 cm	21 - 31 cm	26 - 37 cm
Seat width	13 - 22 cm	19 - 28 cm	27 - 37 cm
Lower leg length	13 - 22 cm	19 - 29 cm	19 - 35 cm
Back rest height (seat to top)	26 - 34 cm	33 - 41 cm	40 - 53 cm

### MANUFACTURER

Timoteos Oy  
Palorannantie 40  
28660 Pori, FINLAND

VAT: FI10408720  
email: timoteos@timoteos.fi

### INTERNATIONAL SALES

Mikko Pirkola  
tel. +358 50 385 6203  
email: sales@timoteos.fi